Step-by-step instructions

1.0 Importing CAD libraries for prosthetic components into the exocad® CAD software

Note: CAD libraries for implant analogs see page 4

Previous versions of existing CAD libraries for prosthetic components (incl. from third party providers) should be deleted before importing the new library in order to minimize data errors.

Removing previous versions of the library

1. Choose the local hard disk (C:) within the Windows Explorer

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3. Choose folder *exocad® DentalCAD* for the currently used version of the CAD

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7. Choose and delete folder with the currently used CAD library (in example *CAMLOG_Titanium_base_CADCAM_for_exocad*)

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After removing previous versions of the library or libraries, the import of the current version (of the library) can be started.

1. Copy folder as part of the downloaded zipfile from the CAMLOG/iSy website

 camlog-camlog-titanium-base_CADCAM_bridge_2016-01-25

in folder *implant*.

(positioned in the same path as described for removing the former version of the library – see previous steps # 1 to 6)

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Important:
Do not import a
Zip archive – only
unzipped folders
can be used!

2. Double-check that the proper library has been imported in the *implant* folder

 caastar-n-rs
 cadstar-s-bl
 cadstar-s-syn
 camlog-camlog
 camlog-camlog-titanium-base_CADCAM_bridge_2016-01-25
 camlog-conelog
 ctech-bl
 ctech-el

3. Close folder.


The exocad® CAD software can now be used with the current CAMLOG library

If you have another folder hierarchy please contact your exocad® reseller.

Step-by-step instructions

2.0 Importing CAD Implant analog libraries into the exocad® CAD software


Disclaimer: Modelcreator is available in exocad® software version EB5164 and above. In combination with Implant analogs the following message appears (status March 2017):

To import Implant analog libraries, follow the steps as described in chapter 1.0, but be aware they are saved in a different folder.

Implant analog libraries must be saved at

Exocad > DentalCAD > DentalCadApp > library > modelcreator > implants

If you have another folder hierarchy please contact your exocad® reseller.

Important:

Please select appropriate implant analog (either NT-Trading or CAMLOG). Only one implant provider can be saved at " ... modelcreator > implants "

Select either DIM analogs from NT-Trading...

Computer Name: BSLN038 > OSDisk (C:) > Programme > Exocad > DentalCAD > DentalCADApp > library > modelcreator > implants >

Name	Änderungsdatum	Typ	Größe
 CAMLOG 3.3 one-piece Ti-CAM Blank typ IAC-ME NT DIM Analog	25.04.2017 15:51	Dateiordner	
 CAMLOG 3.8 one-piece Ti-CAM Blank typ IAC-ME NT DIM Analog	25.04.2017 15:51	Dateiordner	
 CAMLOG 4.3 one-piece Ti-CAM Blank typ IAC-ME NT DIM Analog	25.04.2017 15:51	Dateiordner	
 CAMLOG 5.0 one-piece Ti-CAM Blank typ IAC-ME NT DIM Analog	25.04.2017 15:51	Dateiordner	

...or implant analogs from CAMLOG for CAMLOG and CONELOG

Computer Name: BSLN038 > OSDisk (C:) > Programme > Exocad > DentalCAD > DentalCADApp > library > modelcreator > implants

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Please select your choice of analog(s) from the following libraries:

CAMLOG_CAM-Ti-Blank_type_IAC_ME_CAD-Library_exocad_DE

CAMLOG_Ti-Base_CAD-Library_exocad_DE

CONELOG_CAM-Ti-Blank_type_IAC_ME_CAD-Library_exocad_DE

CONELOG_Ti-Base_CAD-Library_exocad_DE